



# REAL-TIME RT-PCR VS TRADITIONAL METHODS IN DRINKING WATER – A EUROPEAN STUDY

J. BRANDÃO\*

IN REPRESENTATION OF LIVIA GOMEZ, JOÃO BRANDÃO, ANNA NAVARRO, ANTONINO MARTINES, ISABELLA SANSEVERINO, ELENA PORCEL-RODRIGUEZ, DIMITAR MARINOV, OKSANA ABLEITNER, FRANZ ALLERBERGER, ADRIE AT SMA, ROSSELLA BRIANCESCO, ARIJANA CENOV, MARIANNA CÍCHOVÁ, ANNA MARIA COCCIA, SIMONETTA DELLA LIBERA, KATRIEN DE MAEYER, HENDRIK JAN DE VRIES, GOFFE ELSINGA, LEO HEIJNEN, ANNA-MARIA HOKAJÄRVI, VERONIKA JANSKÁ, JERKO HRABAR, MONIKA HUSE-KUTOWSKA, SLAVEN JOZIĆ, STEFANIE KLATOVSKY, CHRISTIAN KORNSCHÖBER, KATHARINA KOHLS, ADRIE KROM, KATRIN LUDEN, CLAUDIA MIKULA-PRATSCHKE, TARJA PITKÄNEN, MATEJA POJE, MILOSLAVA PROKŠOVÁ, MATJAZ RETELJ, BEATA ROZWADOWSKA, ANNASTIINA RYTKÖNEN, ANNIINA SAREKOSKI, SABINE SCHLAGER, BURKHARD SPRINGER, IVA STOJAN, ELINE STROOBACH, NEVEN SUČIĆ, ELISABETE VALÉRIO, GABI VAN PELT, JEROM VRANKEN, ANA VRDOLJAK TOMAŠ, DARIJA VUKIĆ LUŠIĆ, EWA WŁOS, GERHARD WUBBELS, MARIA ZIOMBSKA, KATJA ZELENIK, MARCIN ZIUZIAKOWSKI, AND TERESA LETTIERI

\*NATIONAL INSTITUTE OF HEALTH DR. RICARDO JORGE AND CE3C – CENTER FOR ECOLOGY, EVOLUTION AND ENVIRONMENTAL CHANGES & CHANGE - GLOBAL CHANGE AND SUSTAINABILITY INSTITUTE

<https://orcid.org/0000-0001-7553-227X>  
Email: [joao.brandao@insa.min-saude.pt](mailto:joao.brandao@insa.min-saude.pt)



# STUDY PUBLICATION

Environmental Research 279 (2025) 121786



ELSEVIER

Contents lists available at ScienceDirect

Environmental Research

journal homepage: [www.elsevier.com/locate/envres](http://www.elsevier.com/locate/envres)



## Application of a real-time reverse transcription polymerase chain reaction for rapid detection of *Escherichia coli* in drinking water: an EU representative study

Livia Gómez <sup>a, b</sup>, João Brandão <sup>b, c</sup>, Anna Navarro <sup>c</sup>, Antonino Martinez <sup>d</sup>, Isabella Sanseverino <sup>a</sup>, Elena Porcel-Rodríguez <sup>a, b</sup>, Dimitar Marinov <sup>e</sup>, Oksana Ableitner <sup>f</sup>, Franz Allerberger <sup>f</sup>, Adrie Atsma <sup>g</sup>, Rossella Briancesco <sup>b, h</sup>, Arijana Cenov <sup>i</sup>, Marianna Cíchová <sup>j</sup>, Anna Maria Coccia <sup>h</sup>, Simonetta Della Libera <sup>h</sup>, Katrien De Maeyer <sup>k</sup>, Hendrik Jan de Vries <sup>l</sup>, Goffe Elsinga <sup>m</sup>, Leo Heijnen <sup>m</sup>, Anna-Maria Hokajärvi <sup>n, o</sup>, Veronika Janská <sup>j</sup>, Jerko Hrabar <sup>o, p, b</sup>, Monika Huse-Kutowska <sup>q</sup>, Slaven Jozić <sup>o</sup>, Stefanie Klatovsky <sup>f</sup>, Christian Kornschöber <sup>f, b</sup>, Katharina Kohls <sup>r</sup>, Adrie Krom <sup>g</sup>, Katrin Luden <sup>r</sup>, Claudia Mikula-Pratschke <sup>f</sup>, Tarja Pitkänen <sup>n, t, b</sup>, Mateja Poje <sup>u</sup>, Miloslava Prokšová <sup>j</sup>, Matjaz Retelj <sup>v</sup>, Beata Rozwadowska <sup>w</sup>, Annastiina Rytönen <sup>n, t</sup>, Anniina Sarekoski <sup>n, t, b</sup>, Sabine Schlager <sup>f, b</sup>, Burkhard Springer <sup>f</sup>, Iva Stojan <sup>o</sup>, Eline Stroobach <sup>g</sup>, Neven Sučić <sup>i, b</sup>, Elisabete Valério <sup>b</sup>, Gabi van Pelt <sup>x</sup>, Jerom Vranken <sup>y, b</sup>, Ana Vrdoljak Tomáš <sup>o, b</sup>, Darija Vukić Lušić <sup>i</sup>, Ewa Włos <sup>z</sup>, Gerhard Wubbels <sup>l</sup>, Maria Ziomska <sup>z</sup>, Katja Zelenik <sup>v, b</sup>, Marcin Ziuziakowski <sup>q</sup>, Teresa Lettieri <sup>a, \* b</sup>

<sup>a</sup> European Commission, Joint Research Centre (JRC), Ispra, Italy

<sup>b</sup> INSA National Institute of Health Doutor Ricardo Jorge, CE3C - Centre for Ecology, Evolution and Environmental Changes, Faculdade de Ciências da Universidade de Lisboa, Campo Grande 016, 1749-016, Lisboa, Portugal

<sup>c</sup> ZeClinics Carrer de Laureà Miró, 408-410, Barcelona, Spain

<sup>d</sup> Centro Ricerche Gruppo CAP, Segrate, Italy

<sup>e</sup> ARHS Developments S.A., 4370, Belvaux, Luxembourg

<sup>f</sup> AGES Agency for Health and Food Safety, Graz, Austria

<sup>g</sup> Vitens Water Expertise Centre Leeuwarden, the Netherlands

<sup>h</sup> ISS Italian National Institute of Health, National Center for Water Safety (CeNSia) Rome, Italy

<sup>i</sup> Teaching Institute of Public Health of Primorje Gorski Kotar County, Rijeka, Croatia

<sup>j</sup> WRI Water Research Institute, Bratislava, Slovakia

<sup>k</sup> Pidpa, Antwerp, Belgium

<sup>l</sup> WLN, Glimmen, the Netherlands

<sup>m</sup> KWR Water Research Institute, Nieuwegein, the Netherlands

# CONTEXT & MOTIVATION

- EU Drinking Water Directive 2020/2184 requires zero *E. coli* per 100 mL
- Traditional detection: culture-based methods (21–24h)
- Need for rapid detection in case of contamination events, post-natural disasters scenarios and massive social gatherings

# STUDY OBJECTIVE

- Validate a real-time RT-PCR method for *E. coli* detection on a pan-european scale
- Organised by the Joint Research Centre (JRC) of the European Commission
- Involved 19 labs across 11 EU Member-States

# PARTICIPATING LABORATORIES

11 countries+Eur  
Comission (19  
laboratories): AT,  
BE, HR, DE, FI, IT,  
NL, PL, PT, SI, SK,  
EC

Laboratories participating in the interlaboratory studies. Member States (MS), organisation, laboratory and methods used for RNA extraction and real-time reverse transcription PCR (RT-PCR) are included. N/A, not applicable.


Member State	Organisation	Laboratory	Extraction method	Real-time RT-PCR
Austria (AT) <sup>a</sup>	AGES Austrian Agency for Health and Food Safety	Reference Laboratory for <i>E. coli</i> including VTEC	ZymoBIOMICS DNA/RNA Mini kit	SensiMIX II Probe kit Roche, LightCycler 480
Belgium (BE)	Pidpa	N/A	NucLisSens extraction Biomerieux (no bead beating)	One-step mix - Evocscript Roche Biorad CFX96
	De Watergroep, Heverlee <sup>2</sup>	N/A	Qiagen RNeasy Powerwater (50) Lot 172037190	BioRad CFX96 Deep Well Real-Time System + C1000 Touch Thermal Cycler
Croatia (HR)	IPH Institute for Public Health in the Primorje-Gorski Kotar County	N/A	ZymoBIOMICS DNA/RNA Mini kit	Meridian bioscience Appliedbiosystems Quantstudio 5
	IOP Institute of Oceanography and Fisheries	Laboratory of Microbiology	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Roche, LightCycler 480
Germany(DE)	NLGA Public Health agency of Lower Saxony	N/A	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Bio-RAD CFX96
Finland (FI)	Finnish Institute for Health and Welfare	Laboratory of Water Microbiology, Microbiology Unit	ZymoBIOMICS DNA/RNA Miniprep Kit Modifications to the RNA extraction protocol. Total nucleic acid extraction, followed by TURBO DNA-free™ Kit to remove DNA to obtain pure RNA.	SensiMIX II Probe kit QuantStudio™ 6 Flex Real-Time PCR System, 96-well (Applied Biosystems)
Italy (IT)	Italian National Institute of Health	National Center For Water Safety	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Bio-RAD CFX96 C1000 Touch
The Netherlands (NL)	WLN	N/A	ZymoBIOMICS DNA/RNA Miniprep Kit	HawkZ05 Fast One-step RT-PCR Kit Quantstudio 5 Thermo Fisher Scientific
	Vitens N.V. Water Expertise Centre	N/A	NucLisSens extraction Biomerieux	SensiMIX II Probe kit CFX96 Touch Real-Time PCR
	AQZ	N/A	NucLisSens extraction Biomerieux	SensiMIX II Probe kit QuantStudio™
	Het Waterlaboratorium	N/A	NucLisSens extraction Biomerieux Filters and RNA shield fluid were transferred into Lysisbuffer after which our own protocol was followed.	Roche EvoScript RNA Probes Master Bio-Rad CFX96 C1000 Touch
	KWR Water Research Institute	N/A	NucLisSens extraction Biomerieux	Roche EvoScript RNA Probes Master Bio-Rad CFX Opus 96 Real-Time PCR System
Poland (PL)	Voivodship Sanitary and Epidemiological Station, Unit of the State Sanitary Inspection	Interdisciplinary Laboratory of Molecular Diagnostics in Katowice	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit LightCycler 480 II
		Laboratory of Microbiological Analysis in Bydgoszcz	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Rotor-Gene® Q MDx CE
Portugal (PT)	National Institute of Health Doutor Ricardo Jorge	Laboratory of Environmental and Food Analysis in Olsztyn	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Roche LightCycler96
		Department of Environmental Health	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Bio-Rad CFX96
Slovakia (SK)	WRI Water Research Institute	National Water Reference Laboratory in Slovakia	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Brilliant III Ultra-Fast QPCR Master Mix
Slovenia (SI)	NLZOH National Laboratory of Health, Environment and Food	N/A	ZymoBIOMICS DNA/RNA Miniprep Kit	AriaMX, Agilent, real-time PCR SensiMIX II Probe kit Bio-Rad CFX96
European Commission (EC) <sup>b</sup>	Joint Research Centre (JRC)	Laboratory of Molecular Ecology, Water and Marine Resources	ZymoBIOMICS DNA/RNA Miniprep Kit	SensiMIX II Probe kit Applied Biosystems 7900 HT

<sup>a</sup> Participated only in the first interlaboratory study.

<sup>b</sup> Participated only in the second interlaboratory study.


# METHODS OVERVIEW (1<sup>ST</sup> STUDY -2021)

- The JRC prepared and shipped artificially contaminated water samples (0–10 CFU/100 mL).
- The participating laboratories analysed the samples upon arrival using both the reference and alternative methods




# RESULTS (1<sup>ST</sup> STUDY)

Boxplot of the obtained Ct values by participating laboratories for each contamination level (0, 1, 3, 5 and 10 CFU/100 mL)



Circles indicate outliers according to the R programming language for Statistical Computing. Outliers are defined by the interquartile interval (Q3 - Q1) as values outside the range  $[Q1 - k(Q3 - Q1), Q3 + k(Q3 - Q1)]$  where  $k=1.5$ ,  $Q1=25$ th percentile and  $Q3=75$ th percentile of data

# RESULTS (1<sup>ST</sup> STUDY): Z-SCORES OF THE ALTERNATIVE METHOD




## KEY RESULTS (1<sup>ST</sup> STUDY)

- 18 data sets produced by the participating laboratories.
- 4/18 laboratories reported problems during RNA extraction step and were excluded from the analysis.
- Sensitivity analysis showed lower sensitivity (65.4%) compared to the reference method (90.4%).
- Participating laboratories using polycarbonate filters instead of mixed cellulose filters achieved higher sensitivity of the alternative method (89.2%) compared to the reference method (87.7%).

# METHODS OVERVIEW (2<sup>ND</sup> STUDY -2023)


Artificially contaminated water samples (0–10 CFU/100 mL) were filtered and sent out to the participating laboratories

- The JRC prepared the samples and carried out all filtrations.
- A set of already filtered samples was sent to the participating laboratories that analysed them by the alternative method.
- The JRC analysed all samples by the reference method.



## RESULTS (2<sup>ND</sup> STUDY)


Variability of the water sample analysis by the reference method:  
boxplot CFU/ 100 mL by contamination level (0, 1, 3, 5 and 10 CFU/100 mL)



Circles indicate outliers according to the R programming language for Statistical Computing. Outliers are defined by the interquartile interval (Q3 - Q1) as values outside the range  $[Q1 - k(Q3 - Q1), Q3 + k(Q3 - Q1)]$  where  $k=1.5$ ,  $Q1=25$ th percentile and  $Q3=75$ th percentile of data


## RESULTS (2<sup>ND</sup> STUDY)

Boxplot of cycle threshold (Ct) values by participating laboratories for each contamination level (0, 1, 3, 5 and 10 CFU/100 mL)



Circles indicate outliers according to the R programming language for Statistical Computing. Outliers are defined by the interquartile interval (Q3 - Q1) as values outside the range  $[Q1 - k(Q3 - Q1), Q3 + k(Q3 - Q1)]$  where  $k=1.5$ ,  $Q1=25$ th percentile and  $Q3=75$ th percentile of data


# RESULTS (2<sup>ND</sup> STUDY): Z-SCORES OF THE ALTERNATIVE METHOD




## KEY RESULTS (2<sup>ND</sup> STUDY)

- All laboratories showed performance Z-scores for the considered five contamination levels which are below the critical threshold of 2.
- Acceptability Limits (AL) of the relative level of detection (RLOD) is set at 2.5 according to ISO 16140 part 2 and in this study, the RLOD = 1.79
- Sensitivity – Reference: 97.2%, RT-PCR: 91.1%
- Specificity – Reference: 100%, RT-PCR: 97.7%

# DETECTION TIME COMPARISON



# SENSITIVITY COMPARISON



# CHALLENGES & LEARNINGS

- Lab experience impacts performance
- Critical steps: RNA extraction & filter choice
- Polycarbonate filters perform better than cellulose for nucleic acid extraction

# DISCUSSION – TECHNICAL LESSONS

- Filter type crucial – polycarbonate better for RNA recovery
- Proper lysis (vortex/bead beating) improves sensitivity
- PCR master mix must match equipment (ROX concentration)
- Requires operator training and good lab practice

# DISCUSSION – OPERATIONAL INSIGHTS

- RT-PCR requires training & practice
- Risk of contamination during RNA work
- Standardisation & SOP refinements are key to success

# CONCLUSIONS

- RT-PCR enables rapid detection of *E. coli* in drinking water
- Valid alternative when quick response is needed
- Could complement culture methods, not replace them
- EU-wide study supports adoption with proper validation

# FOR OPEN DISCUSSION

In what context does this molecular method should have an application?

1. as a routine method? (cost/specialised staff)
2. in case of emergency (outbreak/major disaster)?
3. in need of speed?
4. In case of massive social gatherings?
5. some or all of the above?

THANK YOU

